

O₂ Therapy
And
Accessories

Oxygen Therapy



tyco

Healthcare

Oxygen Line

Complete Product Portfolio for Aerosol and Oxygen Therapy

Mask for Aerosol Therapy DAR


Adult

REF. 220/5000

Paediatric


REF. 210/5012

Packaging: 25 pcs.



Mask for Aerosol Therapy Salter Labs

- Mask for Aerosol therapy in clear soft plastic
- Anatomically shaped
- Positioning with Over-the-Ear Style or elastic strap



REF.	DESCRIPTION	PCS. PER PACKAGING
13100	Mask for Aerosol therapy for Adults with fixation - Over-the-Ear Style	50
13101	Mask for Aerosol therapy for Adults with fixation - elastic strap	50
13102	Mask for Aerosol therapy for Paediatrics with fixation - elastic strap	50

Mask for Aerosol Therapy with Nebulizer

Complete with 2 m crush-resistant tubing


Adult

REF. 223/5008

Paediatric

REF. 213/5020

Packaging: 25 pcs.



Oxineb

Complete Nebulizing System from KENDALL range



- Particle size in range of 1 to 4 μm
- Maximum Volume: 8 ml
- Suggested Minimum Volume: 4 ml
- Overspill safety system: 4 ml
- Easy to clean and Setup
- With Oxygen tubing 200 cm

REF.	DESCRIPTION	PCS. PER PACKAGING
21001	OXINEB-N Nebulizer with Adult mask and 200 cm oxygen tubing	50
21002	OXINEB-N Nebulizer with Paediatric mask and 200 cm oxygen tubing	50
21003	OXINEB-N Nebulizer with Mouthpiece and 200 cm tubing	50

Oxygen Mask

Complete with 2m crush-resistant tubing


Adult

REF. 221/5001


Paediatric

REF. 211/5013

Packaging: 25 pcs.



Tracheostomy Mask for Aerosol Therapy



REF. 206/5568


Packaging: 25 pcs.

Venturi Mask


Complete with 2m crush-resistant tubing and Venturi adapter.

Packaging: 25 pcs.

REF. ADULT	REF PAEDIATRIC	CONC.	FLOW LPM	COLOUR
222/5002	212/5014	24%	2	Red
222/5003	212/5015	28%	4	Yellow
222/5004	212/5016	31%	6	White
222/5005	212/5017	35%	8	Green
222/5006	212/5018	40%	10	Blue
222/5007	212/5019	50%	12	Orange
222/5555	212/5554	60%	15	Pink



Mask with Venturi set



With 15 cm flexible tubing, Venturi adapters 2m oxygen crush-resistant tubing and humidification cup

Packaging: 25 pcs.

REF. ADULT	REF PAEDIATRIC	CONC.	FLOW LPM	COLOUR
224/5009	214/5021	24%	2	Red
224/5009	214/5021	28%	4	Yellow
224/5009	214/5021	31%	6	White
224/5009	214/5021	35%	8	Green
224/5009	214/5021	40%	10	Blue
224/5009	214/5021	50%	12	Orange
224/5009	214/5021	60%	15	Pink

Oximask

Mask for Oxygen therapy KENDALL

- Oxygen therapy Mask with 7 different concentrations with Venturi type adapter with regulating O₂% in values:
Base: 24%, 28% (green connector);
High: 35%, 40% (ivory connector)
- Anatomic mask in soft PVC
- Regulating nose-clip fixation
- Elastic strap in Adult version
- With 2 m oxygen tubing



REF.	DESCRIPTION	PCS PER PACKAGING
20001VM	OXIMASK Oxygen Mask - Adult, 7 concentrations, with 2m oxygen tubing	50
20002VM	OXIMASK-TC Oxygen Mask - Adult, 7 concentrations with 2m oxygen tubing and corrugating tube	50
20030VM	OXIMASK-FIX Oxygen Mask - Adult fixed concentration with 2m oxygen tubing and corrugating tube	50
20003VM	OXIMASK-P Oxygen Mask - Paediatric, 7 concentrations, with 2m oxygen tubing	50
20004VM	OXIMASK-SR Oxygen Mask - Adult, w/o regulation, with 2m oxygen tubing	50
20005VM	OXIMASK-SR Oxygen Mask - Paediatric, w/o regulation, with 2m oxygen tubing	50
20101	OXIMASK-T Tracheostomy Mask, 7 concentrations, with 2m oxygen tubing	50
20102	MIX-O-MASK Tracheostomy Mask, 7 concentrations, with 2m oxygen tubing and corrugating tube	50

Venturi Mask Salter Labs

- Transparent
- Anatomic shape, soft plastic
- Oxygen tubing with crush-resistant star shaped inner wall
- System for regulating O₂ concentration with safety block
- Fixation with elastic strap


REF.	DESCRIPTION	PCS PER PACKAGING
13815	Venturi Mask - Adult	50

Aerosol Face Tent

Complete with 2m crush resistant tubing

REF. 202/5024

Packaging: 25 pcs.



Aerosol Face Tent Salter Labs

- Clear mask in soft plastic
- Anatomically shaped
- Fixation with elastic strap




REF.	DESCRIPTION	PCS PER PACKAGING
13103	Aerosol Face Tent with elastic strap	50

Aerosol Face Tent with Nebulizer

Complete with 2m crush resistant tubing

REF. 203/5026

Packaging: 25 pcs.




Aerosol And Oxygen therapy Face Tent mask

With corrugated 1,20 m tubing and a "T" piece

REF. 205/5025

Packaging: 25 pcs.



Oxygen Mask for High Concentrations

Complete with 2m crush resistant tubing

Adult

REF. 225/5011

Paediatric

REF. 215/5023

With non-rebreathing uni-directional valve


Adult

REF. 225/5010

Paediatric

REF. 215/5022

Packaging: 25 pcs.



Oxygen therapy Mask Salter Labs

- Clear
- Anatomically shaped, in soft plastic
- Tubing with crush resistant inner lumen star shaped
- Fixation with elastic strap



REF.	DESCRIPTION	LENGHT cm	PCS PER PACKAGING
13110	Adult Medium Concentraion with tubing (crush resitant) and universal connections	213	50
13122	Adult High Concentration with tubing, non-rebreathing valve and bag, universal connection	213	50
13113	Paediatric Medium Concentration with tubing with connections	213	50

Medication Nebulizer Set Salter Labs

- Nebulizes 3 ml within 10 minutes in vertical or horizontal position based on Venturi principle
- Easy assembling and cleaning
- Larger surface area provided by unique convex cone design to allow total medication volume usage
- Fixation of the mask with elastic strap




REF.	DESCRIPTION	LENGHT cm	PCS PER PACKAGING
13950	Adult with nebulizer and universal connection tubing	213	50
13951	Adult with nebulizer and thread-grip connection tubing	213	50
13962	Adult with Mask and thread-grip connection tube	213	50
13963	Paediatric with Mask and thread-grip connection tube	213	50

Econeb


Complete with nebulizer, "T" Mouth piece, Bacterial/Viral filter Barrierbac "S" and 2 m crush-resistant tubing

REF. 260/5578

Packaging: 25 pcs.



Nebulizer Set



Complete with Nebulizer, "T" piece and 2 m crush-resistant tubing

REF. 260/5032


Packaging: 25 pcs.

Nebulizer

Nebulizes particles in range 2 - 5 μ m

REF. 260/5027

Packaging: 25 pcs.



Humidification Jar



REF. 253/5030


Packaging: 25 pcs.

Safety Tube

2 m

REF. 251/5028

Packaging: 25 pcs.




Safety Tube Salter Labs

- Crush resistant oxygen tubing with star shaped inner lumen
- Available in different connection versions universal on thread-grip connections
- In clear green colour



REF.	DESCRIPTION	LENGHT cm	COLOUR	PCS PER PACKAGING
13400	Safety tube for oxygen and 2 universal connectors	213 cm	Clear Green	50
13409	Safett tube for oxygen and 2 thread-grip connectors	213 cm	Clear Green	50
13444	Adaptor for coupling oxygen tubing	-	-	50

Nasal Cannulas



Complete with 2 m safety tubing

Prematures

REF. 251/5572

Neonatal

REF. 251/5544

Infant

REF. 251/5571


Adult

REF. 251/5029

Packaging: 25 pcs.

Nasal Cannulas Salter Labs

- Nasal or eyeglass cannulas
- Transparent
- Anatomically shaped, soft plastic
- Crush resistant tubing with star shaped inner lumen



REF.	DESCRIPTION	PCS PER PACKAGING
13300	Nasal cannula - Adult - standard with oxygen tubing 213 cm and universal connections	50
13308	Nasal cannula - Adult - High Flows with oxygen tubing 213 cm thread connections	50
13309	Nasal cannula - Adult - Low flows with oxygen tubing 213 cm thread connections	50
13320	Nasal cannulas - Prematures - with 3 mm nasal prongs distance	50
13321	Nasal cannulas - Prematures - with 4 mm nasal prongs distance	50
13322	Nasal cannulas - Infants - with 4.5 mm nasal prongs distance	50
13323	Nasal cannulas - Paediatric - with 5 mm nasal prongs distance	50
13324	Nasal cannulas - Paediatric - with 7 mm nasal prongs distance	50

Oxigafa

Nasal Cannulas KENDALL

- Nasal or eyeglass cannulas
- Anatomically shaped
- Soft PVC
- Oxygen connecting tubing in 190 cm length




REF.	DESCRIPTION	PCS PER PACKAGING
20202	OXIGAFA Nasal Cannula with 190 cm oxygen tubing	100

Pentaspir

REF. 300/6367

Single packaging, sterile: box of 10 pcs.




Catheters and Oxygen connecting tubes Argyle

- In tubes with universal connection at green bubble it is possible to cut bubble expansion

every 91 cm in order to originate funnel connectors, tapered or integral

REF.	DESCRIPTION	LENGHT cm	P.D.	COLOUR	PCS PER PACKAGING
8888-230201	Oxygen tubing, univ.con. at green bubble	30,5 m	3 mm	Clear Green	1
8888-230409	Oxygen joining tubing with tunnel connector	180 cm	3 mm	Clear Green	50
8888-232017	Oxygen Catheter	220 cm	10 Ch	Clear Green	50
8888-232025	Oxygen Catheter	220 cm	14 Ch	Clear Green	50
8888-232215	Oxygen Catheter	40 cm	10 Ch	Clear Green	50

PE FLEX TUBE



REF. 290/5035

Packaging: 50 m Role.

All products are registered trademarks of Tyco Healthcare. © 2001 Tyco Healthcare. All rights reserved.

ECE-CAT-TI4/E

tyco
Healthcare

Tyco Healthcare Italia SpA
Via Rivoltana, 2/D
20090 SEGRATE (MI)
ITALY
Tel: +39 02 703171
Fax: +39 02 70317317
www.tycohealthcare.it
info_italy@tycohealth.com

Tyco Healthcare Switzerland Ltd.
Roosstrasse 53
Postfach
8832 WOLLERAU
SWITZERLAND
Tel: +41 1 786 50 50
Fax: +41 1 786 50 10
www.tycoint.com

Tyco Healthcare ECE
154 Fareham Road
GOSPORT Hampshire
PO 13 0AS
UK
Tel: +44 (0) 1329 224258
Fax: +44 (0) 1329 224086
www.TycoHealth-ECE.com

I-CAT/TI4