

Convective air warming means active warming therapy that provides effective heat transfer over a large surface area and actively maintains and increases body temperature.

WarmTouch™ Warming Unit and Blankets



tyco
Healthcare

MALLINCKRODT

WarmTouch™



WarmTouch™ 5800

- 3 Temperature Settings:
 - High 42-46°C
 - Medium 36-42°C
 - Low 30-34°C
- Operation at room temperature possible.
- Fast Warm-Up (38°C in 30 sec.).
- Smooth durable housing for easy cleaning.
- Quiet operation.
- Multiple mounting options pole, bed rail, cart.
- Electronic overheating prevention.
- For use with WarmTouch™ Blankets for the prevention or treatment of hypothermia.

Technical Data:

Dimensions: 41 cm x 38 cm x 28 cm
(w x h x d)

Weight: 6.8 kg

Power supply: 220 - 240 VAC, 50-60 Hz, 8 A

Ordering Information

Description	REF	QTY
WarmTouch™ 5800	501-5800	1
WarmTouch™ Cart	502-2900	1
Replacement Filter	502-2200	1
Power Cable	502-2500	1
Nozzle	502-2100	1
Hose Assembly	502-2000	1



WarmTouch™ - Blankets for Adults

WarmTouch™ Blankets are available in five convenient sizes to meet all clinical needs. The soft, flexible fabric allows custom draping for any surgical procedure. All blankets contour closely to the patient without billowing up. The precision die-cut holes provide an even airflow. The special outer layer is soft against the patient's skin. No Latex allergens.



CareQuilt™ Full Body Blanket

- For use in Recovery Rooms, A&E, Intensive Care Units, Burn Units, Obstetrics, Outpatient Surgery.
- For optimum warming, blanket should be in direct contact with the patient's skin.
- Chest access slit provides quick and easy access to chest tubes.
- Access slits facilitate checking of arteries and pulses in foot and lower leg.
- Size: 183 cm x 114 cm
- Weight: 310 g
- Also available in paediatric size.

Ordering Information

Description	REF	QTY
CareQuilt™ Blanket	503-0810	25

CareDrape™ Upper Body Blanket

- For use during abdominal, urogenital, and lower-extremity surgery.
- Covers the upper chest and arms during surgery.
- Dual inlet facilitates optimum blower placement and use with two, one or no arm's supports.
- Generous tape length expands draping options.
- Heat loss through head can be minimized by optional head cover.
- Size: 195 cm x 52 cm
- Weight: 290 g

Ordering Information

Description	REF	QTY
CareDrape™ Upper Body Blanket	503-0820	12
Plastic Head Cover	542-0001	50



CareDrape™ Lower Body Blanket

- For use in chest, cardiac and upper extremity surgery.
- Covers lower extremities and lower abdominal region.
- Generous tape length expands draping options.
- Size: 145 cm x 85 cm
- Weight: 305 g

Ordering Information

Description	REF	QTY
CareDrape™ Lower Body Blanket	503-0830	12



CareQuilt™ Paediatric Blanket

- Suitable for children up to 10 years old.
- For use in Burn Units, Intensive Care Units, Outpatient Surgery, Recovery Room.
- Size: 145 cm x 76 cm
- Weight: 165 g

Ordering Information

Description	REF	QTY
CareQuilt™ Paediatric Blanket	503-0840	12



CareDrape™ Paediatric Blanket

- May be used under or over the paediatric patient to help regulate temperature peri-operatively.
- Unique under body design provides unrestricted access to surgical site.
- Clear plastic cover traps warm air.
- X-ray transparent material will not affect or interfere with x-rays.
- Multiple inlets facilitate optimal blower placement.
- Size: 102 cm x 64 cm
- Weight: 101 g

Ordering Information

Description	REF	QTY
CareDrape™ Paediatric Blanket	503-0850	12



Infant:
Blanket under patient



Infant:
Blanket over patient



Pediatric:
Blanket under patient



Pediatric:
Blanket over Patient

WarmTouch™, CareDrape™, and CareQuilt™ are trademarks of Mallinckrodt Inc.
© 2002 Tyco Healthcare. All rights reserved.

M-WT/GB

tyco
Healthcare

Tyco Healthcare UK Limited
154 Fareham Road
GOSPORT
Hampshire
PO13 0AS
UK
Tel: +44 (0) 1329 224226
Fax: +44 (0) 1329 224334

Tyco Healthcare Nederland BV
Hogeweg 105
5301 LL ZALTBOMMEL
The NETHERLANDS
Tel: +31 418 57 66 00
Fax: +31 418 57 67 93
TycoHealthNed@TycoHealth.com
www.tycohealthcare.nl

Tyco Healthcare Regional Export
154 Fareham Road
GOSPORT
Hampshire
PO13 0AS
UK
Tel: +44 (0) 1329 224258
Fax: +44 (0) 1329 224086

Tyco Healthcare Norden AB
Gårdsvägen 18
SE-169 67 SOLNA
SWEDEN
Tel: +46 (0) 8 585 605 05
Fax: +46 (0) 8 585 605 06

Tyco Healthcare Belgium N.V.-S.A.
Generaal de Wittelaan 9/5
B-2800 MECHELEN
BELGIUM
Tel: +32 (0) 15 29 44 50
Fax: +32 (0) 15 29 44 55
tycohealthcare.belgium@tycohealth.com

www.tycohealthcare.com

MALLINCKRODT