



Diagnostic ECG and Ultrasound Supplies

PHILIPS

Dear Customer:

If you're looking for high-quality medical supplies to optimize the performance of your medical devices, you'll find them right here. Whether your focus is on quality patient care, reliability, or cost containment, you can turn to a single source – Philips Medical Supplies. Explore our complete family of supplies, from electrodes to SpO₂ sensors, for every area of care.

With Philips, you can buy your medical equipment and supplies from a single source, simply and easily. Use this guide to find the Philips supplies you need – or go to our website to purchase supplies: <http://www.philips.com/healthcarestore>

Using this guide

Latex Cautions

Philips Medical Supplies currently distributes over 600 products, of which less than 9% contain natural rubber latex or dry natural rubber that could contact humans. Because these substances may cause allergic reactions, products in this guide that contain them are identified with the following symbol:

LATEX

Caution: Products that are identified with this symbol contain dry natural rubber or natural rubber latex that can come into contact with humans and may cause allergic reactions.

To the best of our knowledge, products not identified with this symbol are not known to contain natural rubber latex or dry natural rubber. The latex symbol indicates **only** which **medical supplies** contain natural rubber latex or dry natural rubber that can come into contact with humans and may cause allergic reactions. If you want to know whether a Philips **medical device** contains natural rubber latex or dry natural rubber, refer to the product labeling.

The FDA requires that products be labeled with human contacting latex content information when appropriate. Always refer to product labeling for accurate information on human contacting dry natural rubber or natural rubber latex.

CE Marking

Products labeled with this symbol conform to the European Union (EU) directives as mandated by EU member state legislation and are available for sale in Europe to countries that accept CE marking with restrictions as may be noted. Medical devices labeled with this marking means the product is in conformity with the European

Medical Directive (MDD): 93/42/EEC Council Directive of 14 June 1993 on the approximation of the laws of the Member States concerning medical devices.

In order for a product to be able to be labeled with the CE marking, product quality and manufacturing quality systems must be documented to comply with certain standards. Compliance with the European directives and regulations ensures a standard level of performance for device function, safety, and health protection for users, patients, and others.

The CE symbol provided with some products in the Philips Medical Supplies Reference Guide and online store is intended only as an indication of availability in EU countries and is not intended as a legal product marking.

Countries that accept CE marked products include Austria, Belgium, Bulgaria, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Liechtenstein, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Romania, Russia, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and the United Kingdom. A CE mark does not necessarily indicate that a product is marketed or distributed in all of the countries listed above.

Non-medical products that comply with other harmonized directives may also bear the CE mark.

Ordering Information

To order the medical supplies in this guide or to obtain more information or additional literature, contact your Philips Medical Supplies representative near you or visit our website: <http://shop.medical.philips.com>

Online Ordering

Online ordering is a simple process:

1. Log in using a secure password.
2. Enter the product number and quantity in the shopping cart.
3. Enter your PO or credit card number at check-out.
4. Your order is immediately transferred into our fulfillment system, where you can track it anytime.

If you are new to the web, we can guide you through the process on the phone in minutes. Many new users tell us that the ordering process is so intuitive that they do not need help, even the first time.

You can search for supplies in the online catalog three easy ways: by application category, Philips instrument, or keyword.

If you currently transmit via MMS/EDI or auto fax, the interactive system provides real-time information and services not available through one-way transmissions.

Online ordering is available to U.S. customers at the following url: <http://www.philips.com/healthcarestore>

Diagnostic ECG

Patient Cables and Leads

Patient Cables	7
Patient Data Cables	10
Replacement Lead Sets	11
Extended Leads	13
Electrode-to-Lead Adapters	15

Diagnostic ECG Electrodes

Disposable ECG Electrodes	17
Limb Plate Electrodes	18
Welsh-Bulb Electrodes	19
Replacement Bulbs for Welsh-Bulb Electrodes	19

Cardiograph Paper

Recording Paper	20
-----------------	----

Accessories

Cardiograph Pouches & Cases	23
Cardiograph Batteries	24

Diagnostic Ultrasound

Ultrasound Imaging

Imaging Supplies	25
------------------	----

Imaging Paper

Recording Paper	26
-----------------	----

Transesophageal

TEE Supplies	27
--------------	----

Diagnostic ECG

Patient Cables and Leads

Patient Cables

Product #	Description	Packaging	Use With:
989803129221 CE	Long complete lead set, IEC (4 limb leads, 6 chest leads, 10 color marked rings, 2 lead separators and IFU)	1 lead set	PageWriter Touch: 860284, PageWriter Trim Cardiographs, 860286, 860287, 860288, 860290, 860297
989803129161	Complete lead set, AAMI (4 limb leads, 6 chest leads, 10 color marked rings, 2 lead separators)	1 lead set	PageWriter Touch: 860284, PageWriter Trim Cardiographs, 860286, 860288, 860290
989803129191 CE	Complete lead set, IEC (4 limb leads, 6 chest leads, 10 color marked rings, 2 lead separators)	1 lead set	PageWriter Touch: 860284, PageWriter Trim Cardiographs, 860286, 860288, 860290
989803136951	Standard patient 10 lead set, AAMI	1 lead set	StressVue System: 860295
989803136961	Standard patient 10 lead set, IEC	1 lead set	StressVue System: 860295
989803137661	Large patient 10 lead set, AAMI	1 lead set	StressVue System: 860295
989803137671	Large patient 10 lead set, IEC	1 lead set	StressVue System: 860295



989803129221



989803129161



989803137661



989803137671

Patient Cables *(continued)*

Product #	Description	Packaging	Use With:
989803151631	Complete lead set, AAMI (2 arm leads, 2 leg leads, 2 x 1 into 3 Trident leads, 2 lead separators, 2 Trident Yoke labels (V1-V3, V4-V6), 4 colored clips, 14 colored rings, 6 brown base rings and IFU)	1 lead set	PageWriter TC Series Cardiographs: 860315, 860310
989803151641 CE	Complete lead set, IEC (2 arm leads, 2 leg leads, 2 x 1 into 3 Trident leads, 2 lead separators, 2 banana post adapters, 2 Trident Yoke labels (C1-C3, C4-C6), 4 colored clips, 14 colored rings, 6 white base rings and IFU.	1 lead set	PageWriter TC Series Cardiographs: 860315, 860310
989803151651	Long complete lead set, AAMI (2 long arm leads, 2 long leg leads, 2 x 1 into 3 Trident leads, 2 lead separators, 2 Trident yoke labels (V1-V3, V4-V6), 4 colored clips, 14 colored rings, 6 brown base rings and IFU.	1 lead set	PageWriter TC30 Cardiographs: 860306
989803151661 CE	Long complete lead set, IEC (2 long arm leads, 2 long leg leads, 2 x 1 into 3 Trident leads, 2 lead separators, 10 banana post adapters, 2 Trident yoke labels (C1-C3, C4-C6), 4 colored clips, 14 colored rings, 6 white base rings and IFU.	1 lead set	PageWriter TC30 Cardiographs: 860306
M1713B	Complete leadset, AAMI (4 limb leads, 6 chest leads, 10 color marked rings, 4 shorting plugs, 2 lead separators)	1 lead set	Cardiographs: M1700A, M1701A, M1702A



989803151631



989803151641



M1713B

Patient Cables *(continued)*


Product #	Description	Packaging	Use With:
M1716B CE	Complete leadset, IEC (4 limb leads, 6 chest leads, 10 color marked rings, 4 shorting plugs, 2 lead separators)	1 lead set	Cardiographs: M1700A, M1701A, M1702A
M3702A	Color coded 10-lead patient cable, AAMI Trunk cable and leadset with 1/8" posts	1 lead set	PageWriter Cardiographs: M1770A, M1771A, M1772A, M2662A
M3702C	Color coded 10-lead combination patient cable, AAMI Trunk cable and leadset with 1/8" posts	1 lead set	PageWriter Cardiographs: M1770A, M1771A, M1772A
M3703A CE	Color coded 10-lead patient cable, IEC Trunk cable and leadset with 4mm banana plugs	1 lead set	PageWriter Cardiographs: M1770A, M1771A, M1772A, M2662A
M3703C CE	Color coded 10-lead combination patient cable, IEC Trunk cable and leadset with 4mm banana plugs	1 lead set	PageWriter Cardiographs: M1770A, M1771A, M1772A



M1716B



M3702A



M3703A

Patient Data Cables

Product #	Description	Packaging	Use With:
989803145401	PageWriter Touch USB patient data cable Length: 6.9 feet (2.1m)	1 cable	PageWriter Touch: 860284
989803158481	PageWriter TC70 USB patient data cable , Class A Length: 6.6 feet (2m)	1 cable	PageWriter TC70: 860315
989803164281	PageWriter TC Series USB patient data cable, Class B Length: 6.6 feet (2m)	1 cable	PageWriter TC70: 860315; PageWriter TC30: 860306; PageWriter TC50: 860310
M1719A	PageWriter XL patient data cable Length: 8 feet (2.4m)	1 cable	PageWriter Cardiographs: M1770A, M1771A, M1772A
M1720A	PageWriter XL patient data cable Length: 30 feet (9m)	1 cable	PageWriter Cardiographs: M1770A, M1771A, M1772A



989803145401



989803158481

Replacement Lead Sets

Product #	Description	Packaging	Use With:
989803129141 CE	Limb lead set, AAMI (2 limb leads, 8 color marked rings, 8 colored base rings, and IFU)	1 lead set	PageWriter Touch: 860284, PageWriter Trim Cardiographs, 860286, 860287, 860288, 860290, 860297
989803129151	Chest lead set, AAMI (3 chest leads, 12 color marked rings, 3 brown rings, and IFU)	1 lead set	PageWriter Touch: 860284, PageWriter Trim Cardiographs, 860286, 860287, 860288, 860290, 860297
989803129171 CE	Limb lead set, IEC (2 limb leads, 8 color marked rings, 8 colored base rings, 2 banana plugs and IFU)	1 lead set	PageWriter Touch: 860284, PageWriter Trim Cardiographs, 860286, 860287, 860288, 860290, 860297
989803129181 CE	Chest lead set, IEC (3 chest leads, 12 color marked rings, 3 brown rings, 3 banana adapters, 1 lead separator and IFU)	1 lead set	PageWriter Touch: 860284, PageWriter Trim Cardiographs, 860286, 860287, 860288, 860290, 860297
989803129201 CE	Long limb lead set, IEC (2 limb leads, 8 color marked rings, 8 colored base rings, 2 banana plugs and IFU)	1 lead set	PageWriter Touch: 860284, PageWriter Trim Cardiographs, 860286, 860287, 860288, 860290, 860297
989803129211 CE	Long chest lead set, IEC (3 chest leads, 12 color marked rings, 3 brown rings, 3 banana adapters, 1 lead separator and IFU)	1 lead set	PageWriter Touch: 860284, PageWriter Trim Cardiographs, 860286, 860287, 860288, 860290, 860297



989803129141



989803129151



989803129171



989803129181



989803129201



989803129211

Replacement Lead Sets *(continued)*

Product #	Description	Packaging	Use With:
989803151671 CE	Chest lead set, AAMI and IEC (One 1 into 3 Trident leads, 12x2 colored rings, 3 brown base rings, 3 white base rings (both AAMI and IEC), 3 banana post adapters, 4 Trident Yoke lead description labels (V1-V3, V4-V6, C1-C3, C4-C6) and IFU)	1 lead set	PageWriter TC70: 860315; PageWriter TC30: 860306; PageWriter TC50: 860310
989803151691 CE	Long chest lead set, AAMI and IEC (One 1 into 3 Trident leads, 3x2 colored rings, 3 brown base rings, 3 white base rings (both AAMI and IEC), 3 banana post adapters, 1 lead separator, 4 Trident yoke lead description labels (V1-V3, V4-V6, C1-C3, C4-C6) and IFU)	1 lead set	PageWriter TC70: 860315; PageWriter TC30: 860306; PageWriter TC50: 860310
989803151711 CE	Limb lead set, AAMI and IEC (2 arm leads, 2 leg leads, 5x2 colored rings and 13x2 colored clips for AAMI and IEC, 4 banana post adapters and IFU)	1 lead set	PageWriter TC70: 860315; PageWriter TC30: 860306; PageWriter TC50: 860310
989803151731 CE	Long limb lead set, AAMI and IEC (2 long arm leads, 2 long leg leads, 13x2 colored rings and 5x2 colored clips for AAMI and IEC, 4 banana post adapters and IFU)	1 lead set	PageWriter TC70: 860315; PageWriter TC30: 860306; PageWriter TC50: 860310
M1711A	Limb lead set, AAMI (2 limb leads, 8 color marked rings, 8 colored base rings, and IFU)	1 lead set	PageWriter XL: M1700A, M1701A, M1702A
M1712B	Chest lead set, AAMI (3 chest leads, 12 color marked rings, 3 brown rings, and IFU)	1 lead set	PageWriter XL: M1700A, M1701A, M1702A
M1714A CE	Limb lead set, IEC (2 limb leads, 8 color marked rings, 8 colored base rings, 2 banana plugs and IFU)	1 lead set	PageWriter XL: M1700A, M1701A, M1702A
M1715B CE	Chest lead set, IEC (3 chest leads, 12 color marked rings, 3 brown rings, 3 banana adapters, 1 lead separator and IFU)	1 lead set	PageWriter XL: M1700A, M1701A, M1702A
M1717B CE	Pediatric & Frank lead set, AAMI and IEC (4 chest leads, color marked rings, 4 banana adapters, and IFU)	1 lead set	PageWriter XL: M1700A, M1701A, M1702A



989803151671



989803151711



989803151731

Extended Leads

Product #	Description	Packaging	Use With:
989803148881	16-lead spare parts kit (2 grey and 2 pink marking clips, 5x2 snap-on color rings, 4 banana post adapters, 4 shorting plugs, 2 lead separators and IFU)	1 kit	PageWriter Touch Cardiograph: 860284
989803148931	Upgrade kit 12-15/16 leads, AAMI (4 chest leads, 4 shorting plugs, 2 lead separators, 4 colored clip rings, 5 color marked rings, brown base rings and IFU)	1 kit	PageWriter Touch Cardiograph: 860284
989803149181 CE	Upgrade kit 12-15/16 leads, IEC (4 chest leads, 4 banana post adapters, 4 shorting plugs, 2 lead separators, 5 colored clip rings, 5 color marked rings, white base rings and IFU)	1 kit	PageWriter Touch Cardiograph: 860284
989803149191 CE	Upgrade kit (long) 12-15/16 leads, IEC (4 long chest leads, 4 banana post adapters, 4 shorting plugs, 2 lead separators, 5 colored clip rings, 5 color marked rings, white base rings and IFU)	1 kit	PageWriter Touch Cardiograph: 860284
989803151701	16-lead spare parts kit (9x2 colored marking clips, 4 banana post adapters, 6 shorting plugs, 2 lead separators 2.54 DIA, 2 lead separators 3.4 DIA and IFU)	1 kit	PageWriter TC70: 860315; PageWriter TC50: 860310



989803148931



989803149181



989803149191

Extended Leads *(continued)*

Product #	Description	Packaging	Use With:
989803151751	Upgrade Kit 12-15/16 leads, AAMI (4 additional lead assembly, 6 shorting plugs, 2 lead separators, 6 colored clips, 6 colored rings, 4 brown base rings and IFU)	1 set	PageWriter TC70: 860315; PageWriter TC50: 860310
989803151761 CE	Upgrade Kit 12-15/16 leads, IEC (4 additional lead assembly, 6 shorting plugs, 2 lead separators, 4 banana post adapters, 6 colored clips, 6 colored rings, 4 white base rings and IFU)	1 set	PageWriter TC70: 860315; PageWriter TC50: 860310
989803151771	Upgrade Kit 12-15/16 long leads, AAMI (4 additional long lead assembly, 6 shorting plugs, lead separators, 6 colored clips, 6 colored rings, 4 brown base rings and IFU)	1 set	PageWriter TC70: 860315; PageWriter TC50: 860310
989803151781 CE	Upgrade Kit 12-15/16 long leads, IEC (4 additional long lead assembly, 6 shorting plugs, 2 lead separators, 4 banana post adapters, 6 colored clips, 6 colored rings, 4 white base rings and IFU)	1 set	PageWriter TC70: 860315; PageWriter TC50: 860310



989803151751



989803151761



989803151771

Electrode-to-Lead Adapters

Product #	Description	Packaging	Use With:
989803129231	Alligator clips, AAMI for tab electrodes	10 adapters	PageWriter Trim: 860286, 860287, 860288, 860290; PageWriter Touch: 860284, PageWriter TC30: 860306; PageWriter TC50: 860310; PageWriter TC70: 860315
989803129241 CE	Alligator clips, IEC for tab electrodes	10 adapters	PageWriter Trim: 860286, 860287, 860288, 860290; PageWriter Touch: 860284, PageWriter TC30: 860306; PageWriter TC50: 860310; PageWriter TC70: 860315
989803166031 CE	Clear tab/snap adapter Enables the connection from 3mm - 4mm post banana plug coming from the cardiograph to tab & snap style electrodes.	10 adapters	PageWriter Trim: 860286, 860287, 860288, 860290; PageWriter Touch: 860284, PageWriter TC30: 860306; PageWriter TC50: 860310; PageWriter TC70: 860315
40431B CE	Alligator clip adapter for 1/8" post (AAMI)	10 adapters	PageWriter Trim: 860286, 860287, 860288, 860290; PageWriter Touch: 860284, PageWriter TC30: 860306; PageWriter TC50: 860310; PageWriter TC70: 860315



989803129231



989803129241



40431B

Electrode-to-Lead Adapters *(continued)*

Product #	Description	Packaging	Use With:
40432B CE	Alligator clip adapter for 4 mm banana plug (IEC)	10 adapters	PageWriter Trim: 860286, 860287, 860288, 860290; PageWriter Touch: 860284, PageWriter TC30: 860306; PageWriter TC50: 860310; PageWriter TC70: 860315
40475A CE	Snap lead, adapter AAMI for conversion from 1/8" post to snap electrode connector	10 adapters	For all 1/8 inch post terminated leads
40498E CE	Snap Lead Adapter, IEC for conversion from 4mm banana to grabber electrode connector	10 adapters	For all 4mm banana plug leads
M1718A	Shorting plug set (4 shorting plugs, 10 banana adapters, 2 lead separators)	1 set	PageWriter XL: M1700A, M1701A, M1702A
M2254A CE	Tab electrode adapter to lead wire. Connects tab style electrodes to cardiograph lead wires. Fits on 1/8" post and 4mm banana plug style leads.	10 adapters	For 1/8 inch post and 4mm banana plug leads



40432B



40475A



40498E



M1718A



M2254A

Diagnostic ECG Electrodes

Disposable ECG Electrodes

Product #	Description	Packaging	Use With:
989803149901 CE	Pediatric tab electrode 1.38" x 9/16" (14mm X 34mm)	10/card, 10 cards/pack, 1000/case	PageWriter Touch: 860284
13941E CE	Adult cloth ECG electrode for long-term and holter monitoring 2-1/8" (54mm) diameter	30/pack = 300/case	DigiTrak XT: 860322
13942E CE	Clear plastic pre-gelled electrode. Ideal for stress testing and monitoring 2-1/8" (54mm) diameter	30/pack = 300/case	StressVue: 860295; DigiTrak XT: 860322
13943B CE	Adult solid gel tab electrode for resting ECG 1-3/8" x 7/8" (35mm x 22mm)	100/pack = 1000/case	PageWriter Cardiographs: 860284, 860286, 860287, 860288, 860290, 860297, 860310, 860315
13944B CE	Wet gel snap foam electrode for resting ECG 1-3/4" x 1-1/8" (45mm x 28mm)	10/pack = 300/case	PageWriter Cardiographs: 860286, 860287, 860288, 860290, 860297, 860310, M1700A, M1701A, M1702A, M2406A
40420A CE	High performance snap, pre-gelled foam electrode for resting ECG 2-1/2" x 1-1/8" (65mm x 28mm)	100/pack = 1000/case	Cardiographs, ECG Systems, 860286, 860287, 860288, 860290, M1700A, M1701A, M1702A, 860284, M2406A
M4606A CE	ECG skin preparation paper	100 preps	n/a



989803149901



13941E



13942E



13944B




40420A



M4606A

Disposable ECG Electrodes (continued)

Product #	Description	Packaging	Use With:
M4612A CE	Solid gel snap Electrode	5/pack = 300/case	DigiTrak XT
M4613A CE	Solid gel snap electrode	30/pack = 300/case	DigiTrak XT



M4612A

Limb Plate Electrodes

Product #	Description	Packaging	Use With:
14030A CE	Rubber limb plate strap for electrode (15")	4 straps	Limb Plate Electrodes: 40424A, 40491E, 40494E
40424A	Limb plate electrode for push-in type connection. (not available in Europe)	4 electrodes	Cardiographs with standard 1/8 inch post end leads
40491E CE	Limb plate electrode for 4mm banana plug leads (push-in connector)	4 electrodes	All cardiographs with standard 4mm banana plug leads
40494E CE	Limb clamp electrode, clothespin type for 1/8" post end leads or 4mm banana plug leads, screw connector.	4 electrodes	All cardiographs with standard 1/8 inch post end or 4mm banana plug leads



14030A



40424A




40491E



40494E

Welsh-Bulb Electrodes

Product #	Description	Packaging	Use With:
40418A (LATEX)	Welsh electrode, 5cc bulb (not available in Europe) Shoulder screw connection. Cup diameter: 15mm	6 electrodes	PageWriter Cardiographs
40421A (LATEX)	Welsh electrode, 5cc bulb (not available in Europe) Push-in connector for 1/8" post cable end. Cup diameter: 15mm	6 electrodes	PageWriter Cardiographs
40490E CE (LATEX)	Welsh electrode, 5cc bulb Push-in connector for 4mm banana plug cable end.	1 electrode	PageWriter Cardiographs



40418A



40490E

Replacement Bulbs for Welsh-Bulb Electrodes

Product #	Description	Packaging	Use With:
40419A (LATEX)	Adult size replacement (5cc) bulb for all 15mm diameter electrodes	12 bulbs	Replacement bulbs for 40418A, 40421A, 40490E
40422A (LATEX)	Small size replacement (2cc) bulb for all 15mm diameter electrodes	12 bulbs	




40419A


Cardiograph Paper

Recording Paper

Product #	Description	Packaging	Use With:
989803137011	StressView thermal paper. Z-fold with red grid and no header. Size: 8.5" x 11"	1 case = 12 packs = 3000 sheets	StressVue: 860295
989803143031 CE	StressView thermal paper. Z-fold with red grid and no header. Size: A4 format	1 case = 12 packs = 3000 sheets	StressVue: 860295
CBD20257	Cambridge Heart Stress Chart Paper	1 case = 10 packs = 3000 sheets	Cambridge Heart: CBD2000A
M1707A CE	PageWriter XL thermal paper. Z-fold with 206mm red grid and 2 7/8" (7.3cm) header. Size: 8.5" x 11"	1 case = 8 packs = 1600 sheets	PageWriter XL Cardiographs: M1700A, M1701A, M1702A
M1708A	PageWriter XL thermal paper. Z-fold with 279mm red grid and no header. Size: A4 format	1 case = 8 packs = 1600 sheets	PageWriter XL Cardiographs: M1700A, M1701A, M1702A
M1709A CE	PageWriter XL thermal paper. Z-fold with 212mm red grid and 6.7cm header. Size: A4 format	1 case = 8 packs = 1600 sheets	PageWriter XL Cardiographs: M1700A, M1701A, M1702A
M1710A CE	PageWriter XL thermal paper. Z-fold with 300mm red grid and no header. Size: A4 format	1 case = 8 packs = 1600 sheets	PageWriter XL Cardiographs: M1700A, M1701A, M1702A



989803137011



989803143031



M1707A



M1708A

Recording Paper *(continued)*


Product #	Description	Packaging	Use With:
M1737A CE	Preprinted with 219mm red grid and 2.4 inch (60mm) header ECG laser printer paper. Size: 8.5" x 11"	1 case = 10 packs = 5000 sheets	TraceMaster View: 860281; ECG Management System: M1730A, M1730B
M1738A CE	Preprinted with 240mm red grid and 2.4 inch (60mm) header ECG laser printer paper. Size: A4 format	1 case = 10 packs = 5000 sheets	TraceMaster View: 860281; ECG Management System: M1730A, M1730B
M2206A CE	PTP brand anti-fade thermal paper for PageWriter XL. Z-fold with 206mm red grid and 2 7/8" (7.3cm) header. Size: 8.5" x 11"	1 case = 8 packs = 1600 sheets	PageWriter XL Cardiographs: M1700A, M1701A, M1702A
M2207A CE	PTP brand anti-fade thermal paper for PageWriter XL. Z-fold with 233mm red grid and 6.7cm header. Size: A4 format	1 case = 8 packs = 1600 sheets	PageWriter XL Cardiographs: M1700A, M1701A, M1702A
M2481A CE	Thermal paper for PageWriter cardiographs. Z-fold with 206mm light red grid and 2 7/8" (7.3cm) header. Size: 8.5" x 11"	1 case = 16 packs = 1600 sheets	PageWriter XL: M1700A, M1701A, M1702A; PageWriter 300/200/100: M1770A, M1771A, M1772A; PageWriter Touch: 860284; PageWriter TC70: 860315; PageWriter TC50: 860310; PageWriter TC30: 860306
M2483A CE	Thermal paper for PageWriter cardiographs. Z-fold with 233mm light red grid and 6.7cm header Size: A4 format	1 case = 16 packs = 1600 sheets	PageWriter XL: M1700A, M1701A, M1702A; PageWriter 300/200/100: M1770A, M1771A, M1772A; PageWriter Touch: 860284; PageWriter TC70: 860315; PageWriter TC50: 860310; PageWriter TC30: 860306



M1737A



M1738A



M2207A

Recording Paper *(continued)*

Product #	Description	Packaging	Use With:
M2485A CE	PTP brand anti-fade thermal paper for PageWriter cardiographs. Z-fold with 206mm light red grid and 2 7/8" header (7.3cm) Size: 8.5 x 11"	1 case = 16 packs = 1600 sheets	PageWriter XL: M1700A, M1701A, M1702A; PageWriter 300/200/100: M1770A, M1771A, M1772A; PageWriter Touch: 860284 PageWriter TC70: 860315; PageWriter TC50: 860310; PageWriter TC30: 860306
M2486A CE	PTP brand anti-fade thermal paper for PageWriter cardiographs. Z-fold with 233mm red grid and 6.7cm header. Size: A4 format	1 case = 16 packs = 1600 sheets	PageWriter XL: M1700A, M1701A, M1702A; PageWriter 300/200/100: M1770A, M1771A, M1772A; PageWriter Touch: 860284; PageWriter TC70: 860315; PageWriter TC50: 860310; PageWriter TC30: 860306
M3707A CE	Thermal paper for PageWriter Trim. Z-fold with light red grid and 2 7/8" (7.3cm) header. Size: 8.5" x 11"	1 case = 16 packs = 1600 sheets	PageWriter Trim: 860286, 860287, 860288, 860290
M3708A CE	Thermal paper for PageWriter Trim. Z-fold with light red grid and 7.3cm header. Size: A4 format	1 case = 16 packs = 1600 sheets	PageWriter Trim: 860286, 860287, 860288, 860290



M3707A



M3708A

Accessories

Cardiograph Pouches & Cases

Product #	Description	Packaging	Use With:
989803136651	DigiTrak Plus pouch	1 pouch	DigiTrak Plus Holter Monitor: 860277
989803137961	StressVue pouch for exchange wireless patient module labeled "X12"	1 pouch	StressVue wireless patient module labeled X12
989803137971	StressVue velcro belt with pouch for wired patient module.	1 pouch	StressVue wired patient module
989803144351	StressVue pouch for exchange wireless patient module labeled "X12+"	1 pouch	StressVue wireless patient module labeled X12+
989803153451	DigiTrak XT pouch	1 pouch	DigiTrak XT Holter Monitor: 860322



989803136651



989803137961



989803137971



989803144351

Cardiograph Batteries

Product #	Description	Packaging	Use With:
989803129131	PageWriter Touch Lithium Ion battery (requires 2 batteries)	1 battery	PageWriter Touch: 860284
989803130151	PageWriter Trim rechargeable lead-acid battery	1 battery	PageWriter Trim Cardiographs: 860286, 860287, 860288, 860290, 860297
989803160981	PageWriter TC series Lithium Ion battery	1 battery	PageWriter TC70: 860315; PageWriter TC50: 860310
989803164611	PageWriter Trim long-life battery	1 battery	PageWriter Trim Cardiographs: 860286, 860288, 860290, 860297
CE M2460A	PageWriter 100/200/200i/300i 6-volt battery	1 battery	PageWriter Cardiographs: M1770A, M1771A, M1772A
M3941A	PageWriter 10/10i Ni-MH battery pack	1 battery	PageWriter 10: M2662A



989803129131



989803130151



989803160981



M3941A

Diagnostic Ultrasound

Ultrasound Imaging

Imaging Supplies

Product #	Description	Packaging	Use With:
13921B	Industrial SVHS Tape, Maxell Super VHS Cassette, 1/2"	10 tapes	77612A, 77622A recorders
40483A CE	Aquasonic gel	12 bottles (8.5 oz/ 250ml)	all ultrasound imaging
40483B CE	Ultrasound gel refill, 5 liter bottle (Refill for 40483A)	1 (5 liter) bottle	all ultrasound imaging



40483B

Imaging Paper

Recording Paper

Product #	Description	Packaging	Use With:
40479A	Plasticized thermal strip chart recorder paper (Kowa) Length: 131' (40m) Width: 6" (15.2cm)	4 rolls	Sonos Ultrasound: M2406A, M2407A
40458A	Plasticized thermal video imaging paper (Mitsubishi) Length: 70' (21m) Width: 8.5" (21.59cm)	4 rolls	77020A, 77020AC, 77065AC, 77065AR, 77570A, 77570B, Ultrasound recorders, 770XXA Opt C70, C71, 77507A Opt C70, 77570A Opt C71, 77570B Opt C70, 77570B Opt C71

Transesophageal

TEE Supplies

Product #	Description	Packaging	Use With:
21110A	TEE Transducer disinfection basin and lid	1 Basin and 1 Lid	Philips TEE Transducers
40487A CE	Transesophageal transducer sheath kit	10 procedures	21362A, 21362B, 21362C, 21363A, 21364A, 21367A, 21369A, 21378A, 21478A, M2406A, TEE Transducers
M1828A CE	Disposable TEE bite guard, adult	24 guards	21362A, 21362B, 21362C, 21363A, 21364A, 21367A, 21369A, 21378A, 21478A, M2406A, TEE Transducers
M2203A CE	MiniTee transducer bite guard	24 guards	Agilent 21365A, M2406A, TEE Transducers, 21381A Mini Multi-Plane TEE, 21366A Pediatric Bi-Plane TEE
M2243A	Tip protector for TEE transducer	24 tip protectors	21364A TEE Transducers
M2273A	Tip protector for TEE transducer	24 tip protectors	21367A, 21369A TEE Transducers
M4604A	Tip protector for Omniplane III	24 Tip Protectors	Omniplane III TEE probes: 21378A, 21478A



40487A



M1828A



M2203A



M2243A

Index

13921B.....	25	989803136961.....	7	M1714A.....	12
13941E.....	17	989803137011.....	20	M1715B.....	12
13942E.....	17	989803137661.....	7	M1716B.....	9
13943B.....	17	989803137671.....	7	M1717B.....	12
13944B.....	17	989803137961.....	23	M1718A.....	16
14030A.....	18	989803137971.....	23	M1719A.....	10
		989803143031.....	20	M1720A.....	10
21110A.....	27	989803144351.....	23	M1737A.....	21
		989803145401.....	10	M1738A.....	21
40418A.....	19	989803148881.....	13	M1828A.....	27
40419A.....	19	989803148931.....	13	M2203A.....	27
40420A.....	17	989803149181.....	13	M2206A.....	21
40421A.....	19	989803149191.....	13	M2207A.....	21
40422A.....	19	989803149901.....	17	M2243A.....	27
40424A.....	18	989803151631.....	8	M2254A.....	16
40431B.....	15	989803151641.....	8	M2273A.....	27
40432B.....	16	989803151651.....	8	M2460A.....	24
40458A.....	26	989803151661.....	8	M2481A.....	21
40475A.....	16	989803151671.....	12	M2483A.....	21
40479A.....	26	989803151691.....	12	M2485A.....	22
40483A.....	25	989803151701.....	13	M2486A.....	22
40483B.....	25	989803151711.....	12	M3702A.....	9
40487A.....	27	989803151731.....	12	M3702C.....	9
40490E.....	19	989803151751.....	14	M3703A.....	9
40491E.....	18	989803151761.....	14	M3703C.....	9
40494E.....	18	989803151771.....	14	M3707A.....	22
40498E.....	16	989803151781.....	14	M3708A.....	22
		989803153451.....	23	M3941A.....	24
989803129131.....	24	989803158481.....	10	M4604A.....	27
989803129141.....	11	989803160981.....	24	M4606A.....	17
989803129151.....	11	989803164281.....	10	M4612A.....	18
989803129161.....	7	989803164611.....	24	M4613A.....	18
989803129171.....	11	989803166031.....	15		
989803129181.....	11				
989803129191.....	7	CBD20257.....	20		
989803129201.....	11				
989803129211.....	11	M1707A.....	20		
989803129221.....	7	M1708A.....	20		
989803129231.....	15	M1709A.....	20		
989803129241.....	15	M1710A.....	20		
989803130151.....	24	M1711A.....	12		
989803136651.....	23	M1712B.....	12		
989803136951.....	7	M1713B.....	8		





**Philips Healthcare is part of
Royal Philips Electronics**

How to reach us

www.philips.com/healthcare

healthcare@philips.com

fax: +31 40 27 64 887

Asia

+852 2821 5888

Europe, Middle East, Africa

+49 7031 463 2254

Latin America

+55 11 2125 0744

North America

+1 425 487 7000

800 285 5585 (toll free, US only)

Please visit www.philips.com/healthcarestore



© 2010 Koninklijke Philips Electronics N.V.
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Printed in The Netherlands.
4522 962 54751 * FEB 2010